



EMERSON EXCHANGE 2025

ACCELERATING INNOVATION



ACCELERATING
INNOVATION

6-1701 Merck's Modern Approach to Modular Pilot OT Integration for Process Automation

Disclaimer

The information and/or opinions expressed in this presentation are those of the authors and do not necessarily represent official policy or permission of Emerson or Emerson Exchange.

Important Reminders


Photography and audio/video recording is not permitted in any session, or in the exhibition areas, without press credentials or written permission from Emerson or Emerson Exchange.

Inquiries should be directed to:
EmersonExchange@Emerson.com



Felix Milman

Sr. Specialist, Engineering at Merck & Co., Inc., Rahway, NJ, USA



Eimantas Macys

Sr. Digital Transformation Engineer, Continua Process Systems



Brian Sauerborn

Digital Transformation Engineering Manager, Continua Process Systems





Continua Process Systems

Continua is a team of 60+ NJ-based OT & Process Experts helping labs, plants, and enterprises successfully digitally transform in the Industry 4.0 era.



Merck & Co., Inc., Rahway, NJ, USA

For more than a century, Merck has been at the forefront of research, bringing forward medicines, vaccines, and innovative health solutions for some of the world's most challenging diseases.

Overview

Felix Milman

Brian Sauerborn

Eimantas Macys

Felix Milman

Goals & Business Need

Solution Overview

Implementation Detail

Results

Goals & Business Needs

Goals & Business Needs: User Requirements

- Digital platform** that promotes **interoperability** between benchtop and control layer
- Flexible platform that integrates both common and **uncommon signal types**
- Scalable infrastructure** that minimizes time to integrate and instance new equipment
- Accessible design that is **easy to install, update, and maintain**

Goals & Business Needs: Operational Flexibility

- An R&D pilot plant requires the ability to adapt process design and process equipment.
- Off-the-shelf control systems may introduce operational, design, and communication constraints:
 - May not be designed for **interoperability** with other control systems, may use legacy communication protocols.
 - Limitations on system design, types of equipment that can be used, amount of external I/O that can be connected.
 - Frequent version upgrades making the existing lab system **outdated and unsupported** by vendor.

Goal: Achieve monitoring and control capability for all pilot plant equipment, while minimizing cable clutter

Common Bioprocess Unit Operations & OEM Systems

Unit
Operation

Cell Culture

ATF

TFF

NFF



Goals & Business Needs: Modular Pilot Plant

- An R&D pilot plant's automation infrastructure should be modular to **manage many processes and equipment**.
 - A pilot plant has several concurrent programs running, each with **different process and equipment requirements**.
 - Based on program requirements, a local automation team needs to **adapt configuration quickly** without hindering plant capacity.
- With diverse sets of equipment and control requirements across the plant, 1:1 connections are difficult to manage and instance.
 - When the process or equipment changes at one station, that change must be **instanced to all stations** in the plant.




Figure. Example of a possible lab signal diagram. In this example, it is difficult to manage existing connections and scale out to new stations.

Goal: Design a solution that allows a local automation team to quickly instance and configure equipment and control across the plant from a single platform.

Goals & Business Needs: Digitizing offline data

- **Offline equipment** is often used in cases where an off-the-shelf system does not fit process or communication requirements.
 - This data is **not integrated into the control system**, cannot be monitored remotely, and is difficult to historize and analyze.
- Offline equipment often constitutes a significant percentage of all pilot and lab scale equipment within a plant.
 - These processes rely on manual data transfer, where **process data is not recorded in or monitored in real-time**.
 - Time-to-market is a critical variable in the pharmaceutical development. **Minimizing data transfer and processing time through digital transformation** is a critical business driver.

Goal: Digitize as much offline lab data as possible, minimizing or eliminating manual data-transfer workflows.

Common Offline Equipment

Equipment



Pumps



Probes



Sensors




Balances





Solution Overview

Solution Overview: Guiding Principles of Development

- Compatibility (with POL and most I/O types)
 - Must be robust enough to accommodate nearly any type of I/O
- Compactness (decrease clutter and noise)
 - Must be compact enough to be tucked away and allow for a clean benchtop, but
 - Must be capable enough to accommodate a high density of I/O
- Expandability/Scalability (adapting to technical development environments in lab spaces)
 - Ease of digitization must be at the forefront of decision making




Solution Overview: Scaled-Out



Implementation Detail



Technology (Hardware)

- Emerson RSTi-EP Series Controllers and Slice I/O
 - Contextualizes I/O into the universal MTP Standard
- Network Address Translator
 - Private-side Network for final control elements
- M12 Bulkheads for all signals, power output, ethernet
 - Industry-standard way for swift, robust, and reliable connection
- Expandability due to design





Technology (Hardware) cont.

- Small-footprint
 - standard size: 9”H, 30”W, 8”D
- Single-wire to Distributed Control System
 - Reduced clutter
 - Mobility
- Reliability and Troubleshooting
 - all mechanical/electrical troubleshooting is done between the equipment and SmartBox




NAMUR/MTP

- NAMUR Data Model Generated via MTP Function Blocks within PLC Logic (Process Equipment Assembly)
 - auto-contextualizes data points into a universal standard
 - Like User Defined Type
- All MTP Objects are exposed via an OPC-UA Server with up to 8 Clients
 - connections to a Process Orchestration Layer (DCS, PCS), historians, and data brokers
- Seamless MTP Import to Process Orchestration Layer
 - auto import utilities for MTP are supported by most Distributed Control Systems



Compatibility (Final Control)

- Traditional I/O
 - Analogue Input/Output
 - Digital Input/Output
 - Relay (Dry/Wet Contact)
- Serial I/O
 - Modbus RTU
 - UART
- TCP I/O
 - Modbus TCP





Results

Results:

A new digital foundation

- **+2500 new 1:1 equipment connections added, including 30,000 I/O points**, significantly expanding automation capabilities.
 - Extensive and unprecedented amount of new 1:1 instrument – DeltaV connections made possible through the Smart IO project.
- Smart IO digitization enables a more **ergonomic** lab space.
 - Smart IO consolidates and eliminates tangled wiring, cabling, and bulky HMI's – decluttering and allowing for efficient use of lab space.
- Create a unified Smart IO platform applicable across scales.
 - Off-the-shelf systems are typically designed for specific batch sizes. Smart IO generalizes the integration so that one format can be applied to both benchtop and large-scale processes in an R&D facility
- Smart IO lays a groundwork for **digital transformation**.
 - Integrating hundreds of offline instruments accomplishes an important step towards digital transformation, allowing us to explore new platforms for automating data transfer, alarming, and batch reporting.

Transformation Snapshot Benchtop-scale Perfusion TFF

Comparison	Previous state	Current State
TFF Control	<ul style="list-style-type: none"> • Recirculation: flow sensor to pump 	<ul style="list-style-type: none"> • Recirculation: flow sensor to pump • Permeate: flow sensor to pump • General: any sensor to any pump
External I/O availability	<ul style="list-style-type: none"> • 4x Analog • 4x Serial RS-232 	<ul style="list-style-type: none"> • +8x Analog • +2x Serial RS-232 • +4x ModBus RS-485 • +3x Ethernet
Equipment types	<ul style="list-style-type: none"> • TFF I/O must be compatible with OEM control systems • Bioreactor I/O limited to RS-232 & analog 	<ul style="list-style-type: none"> • Any equipment to support TFF or Bioreactor operations • Addition of PAT for control and process monitoring
User Experience	<ul style="list-style-type: none"> • In-person TFF and bioreactor control • Remote monitoring via 3rd party software integration 	<ul style="list-style-type: none"> • In-person process control via DeltaV • Remote monitoring and control via DeltaV

Goals & Business Needs: User Requirements



Digital platform that promotes **interoperability** between benchtop and control layer



Flexible platform that integrates both common and **uncommon signal types**



Scalable infrastructure that minimizes time to integrate and instance new equipment




Accessible design that is **easy to install, update, and maintain**

Future State

Future Goal: Self Aware Integration

- With the foundation set, the realization of a future where a box of I/O is plugged in and is plug-and-play capable of being digitized and integrated into any platform is closer than ever
- We aim to achieve a future where this box is plugged into a network and automatically “appear” within OT and IT platforms alike, with the right context.
 - As if the Box was “aware” of its existence with respect to the automation platform living within a network



Find More Information

<https://control-associates.emersonimpactpartners.com/products/control-and-safety-systems/conexion-platform/>

Contacts

Hasan Mohammad | Business Development Manager

hmohammad@control-associates.com

201-398-6837





EMERSON EXCHANGE 2025

ACCELERATING
INNOVATION

Thank You